

**Amendments to the Drawings**

The attached replacement sheets of drawings include changes to original figures 1, 2, 2A and 3. These sheets, which include figures 1, 2A, 2B and 3, replace the original drawing sheets 1, 2, 2A and 3. In figure 1 SEQ ID NOS: have been added and the left-hand legend has been corrected. Figure 2 has been relabeled figure 2A, figure 2A has been relabeled figure 2B, and typographical errors have been corrected in figures 2 and 3.

Attachment: Replacement Drawing Sheets

**Amendments to the Specification:**

Please replace paragraph [0052] with the following:

[0052] One method of detection was hybridization of a specific captured target to oligo coupled microspheres and assaying the complex. The reactions were set up as described below. Any allele to be captured was subject to 2 rounds of ~~Hybridization~~ hybridization. The first round of ~~Hybridization~~ hybridization used different homo and heterozygous DNA and specific oligo coupled bead that recognized a particular sequence. The second round of ~~Hybridization~~ hybridization used another set of beads that recognized a specific sequence within that target which confirmed the presence of the captured allele. However, a single round of hybridization was initially done as a control experiment to test the specificity of the oligo coupled microspheres to different alleles within a target.

Please replace paragraph [0072] with the following:

[0072] A second round of ~~Hybridization~~ hybridization was performed with a second bead set that was complementary to the captured template as a test to confirm the accuracy of the template. The samples were measured on a Luminex 100 flow cytometry instrument after the addition of 120 ng of Streptavidin-Phycoerythrin (SA-PE) to each tube and incubated at the hybridization temperature for another 5 minutes. The amount of fluorescent signal obtained was a true representation of the interaction of the biotin with the SA-PE. This assay was a quantitative one and the amount of positive signal was expressed as the highest number obtained for a given reaction.

Please replace paragraph [0068] with the following:

[0068] Different oligonucleotides for specific polymorphisms of the HLA A Locus were coupled to different bead sets (Luminex) to be used in the hybridization assay. The template that hybridized to the oligo coupled beads was selected to provide perfect sequence homology. Coupling beads to specific oligos was performed according to the manufacturer's instructions (Luminex Corp.). The Luminex bead-probe conjugate were hybridized with PCR

fragments produced above. The sequence of the probes used for separation of allele specific PCR fragments was:

L5'A107A    4AGGTATTTCTACACCTCCGTG  
 L5'A107C    4AGGTATTTCTCCACATCCGTG

Please replace paragraph [0069] with the following:

[0069]            The non-hybridized PCR templates were washed away and the PCR ~~fragment~~ fragments ~~specific~~ specifically hybridized to 5'A107A or 5'A107C were eluted from the Luminex beads. Oligos of different sizes, with and without a spacer (i.e. which contained an additional 20 random bases in the middle of an oligo sequence), were coupled to various bead sets and hybridized to different templates to assay for specificity of different alleles. The numbers in the primer identification correlate to different oligonucleotides coupled to beads and indicate the site of the polymorphism for a specific allele. For example, 107 A or C signifies the site of polymorphism at base 107 where each allele either has an A or a C at position 107.

Please replace paragraph [0073] with the following:

[0073]            The second round of hybridization used other allele-specific Luminex bead-probes as follows:

Luminex bead-probes used to confirm allele specific separation

L5'A107A    4AGGTATTTCTACACCTCCGTG  
 L5'A107C    4AGGTATTTCTCCACATCCGTG  
 L5'A153A    4CTTCATCGCAGTGGGCTAC  
 L5'A153C    4CTTCATCGCCGTGGGCTAC  
 L5'A249T    4GCAGGAGGGTCCGGAGTAT  
 L5'A249G    4GCAGGAGGGGCCGGAGTAT  
 L5'A291C    4GAAGGCCCACTCACAGACT  
 L5'A291G    4GAAGGCCCAGTCACAGACT

Please insert the following new paragraphs after paragraph [0074]

[0074.1]       EXAMPLE 3 – Identification of single nucleotide polymorphic allele information in a highly polymorphic HLA-A genomic DNA molecule using allele specific primers coupled to microspheres

[0074.2]       Genomic HLA-A DNA is amplified by PCR as follows. The amplification will be carried out by mixing following components in 20 µl volume in a PCR microtube:

HLA-A locus specific forward and reverse primers 0.25 µM

(In case single-stranded(ss) PCR fragments need to be prepared, the forward primer can have 4 phosphorothioate bonds at their 5' end)

PCR buffers (MgCL<sub>2</sub>/Amonium Sulfate/Tris)  
Deoxynucleotide Triphosphates mix (dNTP), (100uM)  
Genomic DNA (100 ng)  
Tap DNA Polymerase (1 unit)

[0074.3]       PCR is performed using a PE9600 thermal cycler and following the thermal cycling condition: A mini agarose gel is used to check for successful PCR reaction. Following successful PCR reaction, the PCR primers and dNTPs are degraded according to the teachings of Jingwen Chen *et al.*, 2000 which comprises adding 2 units of Shrimp Alkaline Phosphatase and 4 units of E. Coli exonuclease I to the 20 µl of PCR products and incubating at 37°C for 30 min. The PCR enzymes are inactivated by incubation for 15 min at 99°C.

[0074.4]       Single stranded PCR (ssPCR) products are then prepared as taught by Reynolds, *et al.*, 1997. 12 units of T7 gene 6 exonuclease (Amersham Pharmacia) are added to the 20 µl of PCR products and incubated for 30 minutes at 37°C. The loss of double-stranded (ds) PCR fragments can be determined by agarose gel electrophoresis.

[0074.5]       Heterosequence sites in the ssPCR products are then determined by synthesizing a basic set of 20 allele-specific primers to perfectly match all allelic sequences at 10 most common polymorphic sequence sites of the ssPCR products. The 3' end base position of each primer should correspond to a targeted polymorphic base. The 5' end of primer should

be modified with an amino linker. Each primer is conjugated with a different Luminex MAP microspheres (Fulton, *et al.*, 1997). All 20 allele-specific primers are mixed with 4 $\mu$ l (>50 ng) of ssPCR products in 1x PCR buffer, heated for 2 min at 90°C and annealed for 15 min at 55°C in 50  $\mu$ l volume.

[0074.6] The primer extension reaction is then conducted by adding 1 unit of DNA Taq Polymerase, dATP/Cy3-dCTP/dGTP/dTTP to the mixture and incubating at 65°C for 20 min. The fluorescence signals associated with each kind of microspheres are measured using Luminex 100 Flow Cytometer and all heterosequence sites are determined. Heterosequence sites will be identified, if two allele-specific primers targeted at the same heterosequence site are positive.

[0074.7] The allele specific DNA fragment are then separated as follows. The following components are mixed together in a microtube, with separate microtubes used for each allele:

- Biotinylated allele-specific primer/microsphere conjugate targeted to a heterosequence site (500,000 beads)
- Ligation primer (0.2  $\mu$ M)
- Ligation buffer
- 10 units Taq DNA Ligase
- 20  $\mu$ l of ssPCR products as produced above

[0074.8] The components are heated for 2 min at 90°C to remove primers, followed by 30 min incubation at 37°C. Streptavidin coated beads are added to the ligation mixture, and incubated for 15 min at room temperature. The components are then heated 5 min at 75-85°C and the microspheres spun down in a microcentrifuge. The supernatant will contain the nonligated materials. The temperature should be set at the level that the ssPCR fragment hybridized to the ligated primers will not be separated while the ssPCR fragment hybridized to the non-ligated primers will be separated.

[0074.9] If the PCR product does not match the 3' end sequence of the allele-specific primer, the ligation will not occur. The PCR products will be released from the beads-SNP primer due to their lower  $T_m$ .

[0074.10] If PCR product does match the 3' end sequence of allele-specific primer, the ligation will occur. The ligated SNP-pairing primer will form a strong duplex with PCR product that will have significant higher  $T_m$  to sustain 75°C washing step. The supernatant is then removed and the microspheres washed with 500  $\mu$ l wash buffer (2.5 M tetramethylammonium chloride/0.15% SDS/3mMEDTA/75mM TrisHCL, pH 8.0) by repeating the steps of centrifuging and removing the supernatant. The particular alleles coupled to the microspheres can then be determined by conventional genotyping methods with sequencing being a preferred method.

[0074.11] EXAMPLE 4 – Identification of single nucleotide polymorphic allele information in a highly polymorphic HLA-A genomic DNA molecule using allele specific primers and labeled ddNTPs

[0074.12] HLA-A genomic DNA is amplified as in Example 3. Heterosequence sites are determined as described in Example 3. The allele specific ssPCR products are separated as follows. The following components are mixed together in a PCR tube in 20  $\mu$ l volume:

Allele-specific primer targeted to a heterosequence site (0.2  $\mu$ M)  
Biotin-ddATP (2  $\mu$ M)  
Digoxigenin-ddCTP (2  $\mu$ M)  
4 units Thermal Sequenase  
10  $\mu$ l of ssPCR fragments  
4  $\mu$ l of 5x Thermal Sequenase buffer

[0074.13] For this example, we assume A/C polymorphism at the selected heterosequence site. The amount of ssPCR products can be increased if necessary. Single-base extension reaction is performed by cycling 30-50 times at 94°C for 30 sec, 55° for 30 sec, and 72° for 30 sec. The reaction mix is transferred to a QIAquick Spin Column to remove free

Biotin-ddATP and the primers/PCR fragments are eluted in 30  $\mu$ l Tris buffer in a microcentrifuge tube. 5  $\mu$ l of avidin-conjugated magnetic-beads (Dynabeads M-280 Streptavidin) is added and incubated 15 min at room temperature with occasional mixing. A magnet is applied and the supernatant transferred to a second microcentrifuge tube. Digoxigenin-antibody conjugated magnetic beads are added to the second microcentrifuge tube and incubated 15 min at room temperature with occasional mixing. A magnet is applied and the supernatant is removed. The particular alleles coupled to the magnets can then be determined.

[0074.14]      EXAMPLE 5 – Haplotyping of separated alleles of Examples 3 and 4

[0074.15]      Extension primers for all SNPs are synthesized. The extension primers are designed to be complementary to the single-stranded PCR fragment ending one base short of the polymorphic site (at the 3' end). The 5' end of each primer should be modified with amino linker. For HLA-A haplotyping, a complete set of allele-specific primers could be around 90 different primers. Conjugate each primer with a different Luminex MAP microspheres (Fulton, *et al.*, 1997). 100 different Luminex microspheres could be used to conjugate 100 different primers. The individual microspheres can be identified by its fluorescence color ratio between red and infrared and the fluorescence level. Mix all primer/microsphere conjugates (5,000 microspheres each) in the same microcentrifuge tube containing magnetic beads bound with allele-specific PCR fragments from Examples 3 and 4. Divide equally into 4 PCR tubes. Add 2  $\mu$ M of biotin-ddATP to the first tube, biotin-ddCTP to the second tube, biotin-ddGTP to the third tube and biotin-ddTTP to the fourth tube. Add the rest of the three kinds of ddNTP to each tube. Add 4 units of Thermal Sequenase and 4  $\mu$ l of 5xThermal Sequenase buffer. Bring up the final volume to 20  $\mu$ l. Perform single-base extension reaction by cycling 30-50 times at 94°C for 30 sec, 55°C for 30 sec, and 72°C for 30 sec. Add 10  $\mu$ l of diluted avidin-phycoerythrin solution to each tube. Incubate 10 minutes at room temperature. Measure the green fluorescence signals associated with each primer/microsphere conjugates using Luminex 100 Flow Cytometer. The identity of each beads primer conjugate is determined by the red/infrared fluorescence and the polymorphic nucleotide specificity is determined by the green-fluorescence of phycoerythrin.

The recorded results are haplotype/allele-specific which can be analyzed with a specific genotyping software.

[0074.16]      EXAMPLE 6 - Identification of single nucleotide polymorphic allele information in a highly polymorphic HLA-A genomic DNA molecule using allele specific primers covalently attached to variable weight molecules

[0074.17]      HLA-A genomic DNA is amplified as in Example 3. The allele specific ssPCR products are genotyped as follows. Size tagged SNP extension primers are designed as follows. For each SNP, design a SNP extension primer. The 5' part is consisted of a poly-A tail with a different length of A-bases. The 3' part is 18-22 mer oligonucleotide perfect matching the genomic DNA sequence immediately before the polymorphic base site. The length of A-tail between each SNP primer should differ by 2-4 bases, which allows a clear size distinction by gel-electrophoresis. The first base at the 5' end should be a dye-labeled nucleotide. Each primer by itself could run with slightly different mobility due to different mobility due to different nucleotide composition. This has to be considered when deciding the length of oligo-A-tag added to each primer. 20 individual SNPs with each named as SNP1, SNP2, .....SNP20. 20 SNP extension primers with each named as Pex1, Pex2, .....Pex20. Assume that each SNP extension primer has a similar mobility in a gel electrophoresis. The Pex1 will have 2 adenosine added to its 5' end and Pex2 will have four adenosine added to its 5' end. The last SNP extension primer Pex20 will have 40 adenosine attached to its 5' end. Mix 20 SNP extension primers and verify their clear separation in a ABI 377 sequencer machine. Mix following components together in a microtube in a final volume of 20 µl:

- 20 SNP extension primers with an unique size-tag (Pex1 to Pex20)
- 10 µl of ssPCR products
- ddNTPs with each dideoxynucleotide labeled with different dyes
- 4 units of Thermal Sequenase
- 4 µl of 5xThermal Sequenase buffer

[0074.18]      Perform single-base extension reaction by cycling 30-50 times at 94°C for 30 sec, 55°C for 30 sec, and 72°C for 30 sec. Add 20 µl of 2x gel loading mix. Load 2 µl on



ABI 377. Separate 20 SNP extension primers with single extended base by 8% denaturing PAGE. The nucleotide-specificity at each SNP site is determined by fluorescence signal associated with each SNP extension primer band.